

Introducing artificial intelligence for early Vancomycin Resistant Enterococci culture reads with APAS®.



M. Morales¹, B. DeYoung¹, H. Wisplinghoff², R. Green¹ and S. Giglio¹

¹LBT Innovations, Adelaide, Australia

²Labor Dr Wisplinghoff, Cologne, Germany

BACKGROUND

In clinical diagnostic laboratories, the early detection of target bacteria directly impacts early reporting of clinical results and thus patient care. APAS® technology is an advanced artificial intelligence platform used for automated culture plate interpretation. We monitored the growth of four well characterised VRE strains (*E. faecium* ATCC 700221, *E. faecalis* ATCC 51299, *E. faecium* ATCC 51858, *E. faecalis* NCTC 12201) over a period of 48 hours utilising Thermo Fisher Brilliance VRE Agar and APAS Compact. We aimed to compare how early APAS® could reliably detect VRE positive plates when compared to a routine microbiologist read.

METHOD

Four isolates were used in the study: *E. faecium* ATCC 700221, *E. faecalis* ATCC 51299, *E. faecium* ATCC 51858, and *E. faecalis* NCTC 12201. Each Thermo Fisher Brilliance VRE Agar plate was inoculated with high a concentration of VRE colonies (approximately 1000 CFU/plate) and a low concentration of VRE colonies (approximately 10 CFU/plate). Plates were contained within an APAS® Compact were incubated at 36±1°C for 48 hours and imaged every 5 minutes.

A total of 1,152 images were taken per isolate preparation. APAS® and microbiologist assessment were performed on images at every 30-minute timepoint. The images used for microbiologist assessment were scrambled prior to assessment to remove time bias. The earliest time to detection was defined as the earliest point where APAS® or the microbiologist independently determined an image as positive, and when all the other subsequent images remained positive.

RESULTS

FIGURE 1: VRE Time-lapse video




FIGURE 2: APAS® colony detection at 12 hours






FIGURE 3 : VRE plate images at different time intervals using APAS®



RESULTS

FIGURE 4: Earliest time to detection results



Our results demonstrate that the earliest time that APAS® technology detected consistent VRE growth was approximately 13 hours after incubation and that variability existed in the way microbiologists designated positive or negative plates.

SUMMARY / CONCLUSION

In this study, we demonstrated that a significant reduction in reading time for VRE culture plates (<48 hours as recommended by manufacturer) may be possible when using APAS®.

Our results suggest that APAS® can assist in early reads of VRE culture plates and assist with infection control surveillance.

Future research should focus on demonstrating similar behaviour in wild-type isolates.

As the demand for faster turn around times for critical tests increases, efforts should be made to utilise artificial intelligence technologies for early reads where possible.

CE DISCLOSURES / ACKNOWLEDGEMENTS

The project was internally funded by LBT Innovations.

CONTACT INFORMATION

Dr Steven Giglio PhD, FASM
 Scientific Director, LBT Innovations
 Email: Steven.Giglio@lbtinnovations.com